

Otto Bock®

QUALITY FOR LIFE



B600 Power Wheelchair

For Complex Rehabilitation

The B600 Power Wheelchair for complex rehabilitation is built with quality and function from the ground up. We started with spring mounted wheels and a 6.2 mph rear wheel drive for indoor space maneuverability, then we added a stable, telescoping frame that can be adjusted to the user's center of balance.

It's also loaded with options, and more options mean a better fit. A power seat tilt mechanism offers relief of up to 30°, and a seat elevator lifts the user up to 16 inches for those hard-to-reach cupboards! The B600: ***Elevate your expectations!***



B600 Power Wheelchair



Features and Benefits

- Eye-catching design hides the “guts” of the base for a sleeker look
- Spring mounted front* and rear wheels ensures that the chair absorbs the jolts of driving, not the user
- Telescoping frame adjusts the center of gravity for increased stability
- Rear wheel drive allows the user's to easily navigate around indoor obstacles
- Pull-out battery tray and direct drive make it easy to service

*Optional



Technical Data	
Max. Speed	6.2 mph
Battery	group 24 flat top, gel
Distance Range	approx. 22 miles
Charging Time	approx. 12 hours
Max. Occupant Capacity	264 lbs / 396 lbs*
Product Weight (empty)	253 lbs (with batteries)
Obstacle Climbing Height	60 mm
Product Width	25.4"
Product Length	42.5" (with footrests) 33" (without footrests)
Seat to Floor Height	17.5"-23.5" (in 2" increments)
Seat Width	15"-21.5" (standard seat, infinitely adjustable) 20" (Recaro seat, non-adjustable)
Seat Depth	15"-19.5" (standard seat, incrementally adjustable) 19"-21" (Recaro seat)
Back Angle Adjustment	-9° to 30° (standard seat**) -10° to 30° (Recaro seat)
Armrest Height	9.5"-14"
Turning Radius	32.7" (pivot turn) / 44" (corridor turn)

*with heavy duty option

**two adjustment ranges (-9°/1°/11°/21° and 0°/10°/20°/30°)