

# Reciprocating Gait Orthosis (RGO)

USA 800.328.4058

Fax 800.962.2549

Effective 5/01/2007

## Order Form


Date _____	Account No. _____
Account Name _____	Ship to Address _____
Bill to Address _____	City/State/Zip _____
City/State/Zip _____	Phone No. _____
Fax No. _____	Email Address _____
_____	Authorized Signature _____

Purchase Order No. \_\_\_\_\_

The RGO Hip Joint System is suitable for patients up to 143 lbs. The System is available in three different sizes. 17H100=2 for a body weight of up to 75 lbs. 17H100=1/0 a body weight of up to 143 lbs. For the pelvic widths of 200 mm to 380 mm, pelvic tubes in gradations of 1 cm are available. Important: Measure the clear width between the two trochanters while sitting and in prone position. When measuring, take into account the padding as well as the thickness (X) of the pelvic module to be fabricated.

### Please check the desired combination

Combination	17H100=2 up to 75 lbs	17H100=1 up to 143 lbs	17H100 up to 143 lbs
17Y140=200			
17Y140=210			
17Y140=220			
17Y140=230			
17Y140=240			
17Y140=250			
17Y140=260			
17Y130=270			
17Y130=280			
17Y130=290			
17Y130=300			
17Y130=310			
17Y130=320			
17Y130=330			
17Y130=340			
17Y130=350			
17Y130=360			
17Y130=370			
17Y130=380			



Weight \_\_\_\_\_ lbs

Height \_\_\_\_\_ in

Diagnosis \_\_\_\_\_

Clarify thigh bar offset:  0 mm  5 mm  10 mm (17Y131/143)

14630 28th Avenue North, Minneapolis, MN 55447  
For more information please visit our web site at [www.ottobockus.com](http://www.ottobockus.com)

All specifications and pricing are subject to change without notice.  
©2007 Otto Bock HealthCare LP • 07041285.2 • 5/07

**Otto Bock®**

QUALITY FOR LIFE