

# SERVICE PARTS, ACCESSORIES

## Service Parts: 1048-1079 Additional Pads for Dosi Supports

### Plastic Pad in Fabric Pocket



Art. no.	Suitable for Supports
1030	1064, 1069, 1070, 1079

## Service Parts: 50R20 Dorso Osteo Care

### Back Pad



Art. no.	Size
29R249	S
	M

### Abdominal Support



Art. no.	Size
29R261	universal

### Shoulder Straps



Art. no.	Size
29R260	universal

### Waist Straps



Art. no.	Size
29R262	universal














### Clip Rivets




Art. no.	Contents	Colour
504K8	12-pack	white

# SERVICE PARTS, ACCESSORIES




## Service Parts: 4024 - 4146

Splints									
	<table border="1"> <thead> <tr> <th>Art. no.</th> <th>Plastic splints for 4024, 4025, 4055, 4143, 4146</th> </tr> </thead> <tbody> <tr> <td>4002</td> <td> = soft</td> </tr> <tr> <td>4000</td> <td> = medium</td> </tr> <tr> <td>4001</td> <td> = hard</td> </tr> </tbody> </table>	Art. no.	Plastic splints for 4024, 4025, 4055, 4143, 4146	4002	 = soft	4000	 = medium	4001	 = hard
	Art. no.	Plastic splints for 4024, 4025, 4055, 4143, 4146							
	4002	 = soft							
4000	 = medium								
4001	 = hard								
<table border="1"> <thead> <tr> <th>Art. no.</th> <th>Plastic splints for 4058</th> </tr> </thead> <tbody> <tr> <td>4014</td> <td> = hard</td> </tr> </tbody> </table>	Art. no.	Plastic splints for 4058	4014	 = hard					
Art. no.	Plastic splints for 4058								
4014	 = hard								
<table border="1"> <thead> <tr> <th>Art. no.</th> <th>Aluminum splints for 4024, 4025, 4055, 4143, 4146</th> </tr> </thead> <tbody> <tr> <td>4014</td> <td></td> </tr> </tbody> </table>	Art. no.	Aluminum splints for 4024, 4025, 4055, 4143, 4146	4014						
Art. no.	Aluminum splints for 4024, 4025, 4055, 4143, 4146								
4014									

	<p>The support splints follow the contour of the wrist. They are pre-shaped to provide a functional rest position of the wrist as well as to fit to the shapes of the palm, wrist and forearm. Plastic splints can be custom-adapted with hot water (78 - 80 °C / 167 - 176 °F). The material retains its shape and is environmentally friendly. Available in three hardness ranges. The aluminum splints do not contain heavy metals such as chrome, nickel or zinc. They can be customized to the patient's individual contours.</p>
--	--

## Service Parts: 28P50 Wrist and Thumb Orthoses

Cuffs														
	<table border="1"> <thead> <tr> <th>Art. no.</th> <th>Description</th> <th>Material</th> </tr> </thead> <tbody> <tr> <td>29P47=1</td> <td>Cuff</td> <td>Cotton</td> </tr> <tr> <td>29P47=2</td> <td>Cuff</td> <td>Cotton</td> </tr> <tr> <td>29P47=3</td> <td>Cuff</td> <td>Cotton</td> </tr> </tbody> </table>	Art. no.	Description	Material	29P47=1	Cuff	Cotton	29P47=2	Cuff	Cotton	29P47=3	Cuff	Cotton	
Art. no.	Description	Material												
29P47=1	Cuff	Cotton												
29P47=2	Cuff	Cotton												
29P47=3	Cuff	Cotton												
														
														

## Service Parts: 50K13 Genu Arexa

### Service Set



Art. no.	Side		Size
29K80	left	right	S
	left	right	M
	left	right	L
	left	right	XL
	left	right	XXL

- Consisting of:
- 2 x Upper part pad
  - 2 x Lower part pad
  - 2 x Strap pad
  - 1 x Frame pad
  - 1 x Pad cover
  - 1 x Tibia pad
  - 2 x Textil cover for TechnoGel knee pad
  - 2 x TechnoGel Knee pad
  - 2 x Ergonomic closure
  - 6 x Hook-and-loop Closure Y-design  
Micro hook tape, back-to-back
  - 12 x Plastic pin - upper part
  - 12 x Plastic pin - lower part
  - 2 x Stop set
  - 6 x Strap
  - 4 x Distance ring
  - 2 x Condyle basket
  - 2x Clamp guard

## Service Parts: 28K20/28K21 Genu Arthro

### Hook-and-loop Straps



Art. no.	Strap Width		Length		Colour
	cm	inch	cm	inch	
170Z=400 - 7	3.8	1.5	40	15.8	black
170Z=900 - 7	3.8	1.5	90	35.4	black

### Hook-and-loop Closure Y-Design



Art. no.	Strap Width		Colour
	mm	inch	
29Y161=40 - 7	3.8	1.5	black


# SERVICE PARTS, ACCESSORIES

## Ergonomic Closure



Art. no.	Strap Width	
	cm	inch
170D20	3.8	1.5

## Strap Pads



Art. no.	Strap Width		Length	
	cm	inch	cm	inch
29S29	3.8	1.5	44	17.3

## Frame Pad



Art. no.	Side		Size
	left	right	
29S26	left	right	XS
	left	right	S
	left	right	M
	left	right	L
	left	right	XL

## Foot Stirrup



Art. no.	Side	
	left	right
29K56	left	right

## Measuring Template



Art. no.
743A26