

Date	Account No.
Account Name	Ship to Address
Bill to Address	City/State/Zip
City/State/Zip	Phone No.
Fax No.	Email Address
Physician/Therapist	Authorized Signature

Mark For: Quote Only Please Order

Purchase Order No.: _____

Model Type

Lisa™ Pediatric Postural Mobility System
 Each system includes caster swivel wheel locks, reinforced solid dark blue nylon seat cover, back angle adjustment of 30°, integrated tie-down anchor points, swiveling solid wheels with ball bearings, one-step wheel locks, integrated tip assist, height adjustable foot rest and a 2' standard plastic lap belt in black. Comes standard with Lisa. (E1236) **HCPCS**

<input type="checkbox"/> Lisa – Size 1	HRK3213=01-027-USA	E1236	\$2,075.00
<input type="checkbox"/> Lisa – Size 2	HRK3214=01-027-USA	E1236	\$2,185.00

Seat Insert Option:
 Seat insert comes in "Paint", a special patterned fabric.

Padded Seat Insert w/multi-functional back

<input type="checkbox"/> Lisa – Size 1	HR32135200-334	K0108	\$405.00
<input type="checkbox"/> Lisa – Size 2	HR32145200-334	K0108	\$416.00

Postural Supports

Headrest, padded (blue)

<input type="checkbox"/> Lisa – Size 1	HR32070060-031	E0955	\$255.00
<input type="checkbox"/> Lisa – Size 2	HR32086500-031	E0955	\$270.00

Headrest, for multi-functional back insert (blue)

<input type="checkbox"/> Lisa – Size 1	HR32137700-031	E0955	\$270.00
<input type="checkbox"/> Lisa – Size 2	HR32147700-031	E0955	\$270.00

Armrest (each)

<input type="checkbox"/> Lisa – Size 1 & 2	HR32137100-027	K0019	\$75.00
--	----------------	-------	---------

Lap Belt with Metal Buckle

<input type="checkbox"/> Lisa – Size 1	HR32020030	E0978	\$75.00
<input type="checkbox"/> Lisa – Size 2	HR32030030	E0978	\$86.00

Five Point Harness

<input type="checkbox"/> Lisa – Size 1 & 2	HR32136100	E0960	\$190.00
--	------------	-------	----------

Five Point Harness, for multi-functional back insert

<input type="checkbox"/> Lisa – Size 1 & 2	HR32135400	E0960	\$198.00
--	------------	-------	----------

Postural Supports (continued)

Padded Groin Strap (blue) (with or without seat inserts)

<input type="checkbox"/> Lisa – Size 1	HR32075800-031	K0108	\$148.00
<input type="checkbox"/> Lisa – Size 2	HR32085800-031	K0108	\$152.00

Flip-down Abductor (blue)

<input type="checkbox"/> Lisa – Size 1	HR32137200-031	E0957	\$237.00
<input type="checkbox"/> Lisa – Size 2	HR32147200-031	E0957	\$278.00

Width Reduction Pads, for multi-functional back insert (blue)

<input type="checkbox"/> Lisa – Size 1	HR32135500-031	K0108	\$216.00
<input type="checkbox"/> Lisa – Size 2	HR32145500-031	K0108	\$206.00

Note: The width reduction pads are filled with foam inserts.

Grab Rail (metal, without padding)

<input type="checkbox"/> Lisa – Size 1	HR32137800	K0108	\$88.00
<input type="checkbox"/> Lisa – Size 2	HR32147800	K0108	\$88.00

Grab Rail Padding (blue)

<input type="checkbox"/> Lisa – Size 1	HR32135700-031	K0108	\$53.00
<input type="checkbox"/> Lisa – Size 2	HR32145700-031	K0108	\$52.00

Padded Foot straps (pair)

<input type="checkbox"/> Lisa – Size 1 & 2	HR32020080	E0951	\$40.00
--	------------	-------	---------

Folding Storage Bag

<input type="checkbox"/> Lisa – Size 1	HR32135800	K0108	\$125.00
<input type="checkbox"/> Lisa – Size 2	HR32145800	K0108	\$125.00

Sun/Rain Canopy

<input type="checkbox"/> Lisa – Sizes 1 & 2	HR32136200-094	K0108	\$243.00
---	----------------	-------	----------

Medical Necessities Bag

<input type="checkbox"/> Lisa – Sizes 1 & 2	662T13	K0108	\$42.00
---	--------	-------	---------

14630 28th Avenue North, Minneapolis, MN 55447
 For more information please visit our web site at www.ottobockus.com

All specifications and pricing are subject to change without notice.
 ©2010 Otto Bock HealthCare LP • 05100884.11 • 2/10

Lisa™ Technical Specifications:

Dimensions and Weights	Size 1	Size 2
Load Capacity	110 lbs	132 lbs
Product Weight	30 lbs	31 lbs
Push Handle Height	41"	41"
Overall Width	24"	26"
Min. Overall Width	40"	40"
Max. Overall Width	48"	48"
Seat Depth		
Min. Seat Depth	11.75"	13.75"
Max. Seat Depth	13.75"	15.75"
Seat Depth with multi-functional back	min: 10.25" max: 12.25"	min: 12.25" max: 14.75"
Seat Width	14.2"	16"
Seat Width with multi-functional back	13"	14.5"
Seat Width with multi-functional back with width reduction pads	9.5" / 10.25" / 11"	11" / 11.75" / 12.5"
Seat Height		
Anterior	19.7"	20.5"
Posterior	15"	15"
Seat Angle	22°	22°
Backrest Height	25.6"	27.5"
Seat to Back Angle (4 settings in increments of 7.5°)	90° - 120°	90° - 120°
Footrest Height		
Min. Lower Leg Length	10.5"	10.5"
Max. Lower Leg Length	15.75"	15.75"
Front Wheel Diameter	7"	7"
Back Wheel Diameter	10"	10"
Folding Size	41" x 18" x 17"	41" x 18" x 17"
Turning Radius	55"	57"